DuoSet Technical tip

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PREPARATION & STORAGE

Store unopened kit at 2-8 °C. Do not use past kit expiration date.

REAGENT	PREPARATION	STORAGE OF OPENED/RECONSTITUTED MATERIAL
Capture	Reconstitute with 0.5 mL of PBS	Store at 2-8 °C for up to 8 weeks or aliquot and store at -20 °C to -70 °C in a manual defrost freezer for up to 12 weeks.*
Detection	Reconstitute with 1.0 mL of Reagent Diluent	
Standard	Reconstitute with 0.5 mL of Reagent Diluent	Aliquot and store reconstituted standard at -70 °C for up to 12 weeks.*
Streptavidin-HRP	Dilute with Reagent Diluent	Store undiluted at 2-8 °C for up to 12 weeks. DO NOT FREEZE.*

*Provided this is within the expiration date of the kit.

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Troubleshooting your DuoSet® ELISA *continued*

Problem	Possible Cause	Solution
Poor assay to assay reproducibility	Insufficient washing	<ul style="list-style-type: none"> • See washing procedures • If using an automatic plate washer, check that all ports are clean and free of obstructions
	Variations in incubation temperature	<ul style="list-style-type: none"> • Adhere to recommended incubation temperature • Avoid incubating plates in areas where environmental conditions vary
	Variations in protocol	<ul style="list-style-type: none"> • Adhere to the same protocol from run to run
	Plate sealers reused, resulting in presence of residual HRP which will turn TMB blue	<ul style="list-style-type: none"> • Use fresh plate sealer for each step
	Improper calculation of standard curve dilutions	<ul style="list-style-type: none"> • Check calculations, make new standard curve • Use internal controls
	Buffers contaminated	<ul style="list-style-type: none"> • Make fresh buffers
No signal when a signal is expected, but standard curve looks fine	No cytokine in sample or levels below assay range	<ul style="list-style-type: none"> • Use internal controls • Repeat experiment, reconsider experimental parameters
	Sample matrix is masking detection	<ul style="list-style-type: none"> • Dilute samples at least 1:2 in appropriate diluent, or preferably do a series of dilutions to look at recovery
Samples are reading too high, but standard curve looks fine	Samples contain cytokine levels above assay range	<ul style="list-style-type: none"> • Dilute samples and run again
Very low readings across the plate	Incorrect wavelengths	<ul style="list-style-type: none"> • Check filters/reader
	Insufficient development time	<ul style="list-style-type: none"> • Increase development time
	Coated plates are old and have gone bad	<ul style="list-style-type: none"> • Coat new plates
	Capture antibody did not bind to the plate	<ul style="list-style-type: none"> • Use an ELISA plate (not a tissue culture plate) • Dilute in PBS without additional protein
	Buffer containing FCS used to reconstitute antibodies	<ul style="list-style-type: none"> • Requalify your reagents of choice
Green color develops upon addition of stop solution when using streptavidin-HRP	Reagents not mixed well enough in wells	<ul style="list-style-type: none"> • Tap plate
Edge Effects	Uneven temperatures around work surfaces	<ul style="list-style-type: none"> • Avoid incubating plates in areas where environmental conditions vary
		<ul style="list-style-type: none"> • Use plate sealers
Drift	Interrupted assay set-up	<ul style="list-style-type: none"> • Assay set-up should be continuous – have all standards and samples prepared appropriately before commencement of the assay
	Reagents not at room temperature	<ul style="list-style-type: none"> • Ensure that all reagents are at room temperature before pipetting into the wells unless otherwise instructed in the antibody inserts

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