



# Mesenchymal Stem Cells



# Mesenchymal Stem Cells

Mesenchymal stem cells (MSCs) are defined as multipotent, self-renewing progenitors that can be differentiated into adipocytes, chondrocytes, and osteocytes. Originally identified in mouse bone marrow, MSCs have now been discovered in a variety of species and isolated from numerous tissues including adipose, placental, dental pulp, and umbilical cord. Despite the classical trilineage differentiation that functionally identifies MSCs, these cells have also been shown to differentiate into non-traditional lineages to produce cardiomyocytes, endothelial cells, hepatocytes, and neural cells. To date, the biological properties of MSC identification, differentiation, and function have yet to be confirmed *in vivo*, raising caution for the extrapolation of *in vitro* generated data. Given their rarity, incompletely defined immunophenotype, and localization in multiple organs, studying MSCs *in situ* is not a trivial task. With these challenges in mind, R&D Systems and Tocris Bioscience present tools to assist in the reliable isolation, differentiation, verification, and investigation of MSCs.



# Isolate and Culture

A number of protocols are available to isolate MSCs, however each can vary in the yield, purity, quality, and the ability of the cells to proliferate. It is important to begin experimentation with a verified, homogeneous, multipotent MSC starting population to ensure confidence in differentiation and subsequent data interpretation. With this necessity in mind, R&D Systems offers a variety of products to establish and maintain a homogenous starting MSC population, including purified MSCs, MSC isolation kits, and a range of media for MSC expansion.

| Product  | Description   | Catalog # |
|--|---|-----------|
| Rat Mesenchymal Stem Cells                           | CD90 <sup>+</sup> , CD45 <sup>-</sup> , CD31 <sup>-</sup> | PSC003    |
| MagCelect™ Mouse Mesenchymal Stem Cell Isolation Kit | Yields highly purified MSCs (75%–95%)                     | MAGM212   |
| StemXVivo™ Xeno-Free Human MSC Expansion Media       | Free of non-human animal-derived components               | CCM021    |
| StemXVivo Serum-Free Human MSC Expansion Media       | Serum free expansion media                                | CCM014    |
| StemXVivo Mesenchymal Stem Cell Expansion Media      | Media for human, mouse, and rat MSCs                      | CCM004    |
| CryoDefend™-Stem Cells Media                         | For defined MSC cryopreservation                          | CCM018    |

## Rat Mesenchymal Stem Cells


### Features

- **Verified** – confirmed phenotype of CD90<sup>+</sup>, CD45<sup>-</sup>, CD31<sup>-</sup>
- **Multipotent** – confirmed differentiation to adipocytes, chondrocytes, and osteocytes
- **High purity** – homogeneous population reduces experimental variation
- **Ready to use** – free of mycoplasma and microbial contamination


## MagCelect Mouse Mesenchymal Stem Cell Isolation Kit

### Features

- **Rapid** – isolation within 2 hours
- **Efficient** – high yields (75–95%) of pure MSCs
- **Gentle** – negative selection ensures isolated MSCs are undisturbed
- **Robust** – each kit processes up to 3 x 10<sup>8</sup> cells



**Rat MSCs Show Positive Expression of CD90/Thy1 and Negative Expression of CD45 and CD31.** (A) Rat MSCs (Catalog # PSC003) were stained for the indicated markers after passage 5. CD90/Thy1 was detected using an Anti-Rat CD90/Thy1-PE Antibody, CD45 was detected using an Alexa Fluor® 647-conjugated Anti-Rat CD45 Antibody. (B) CD31 was detected using a Goat Anti-Mouse CD31/PECAM-1 Affinity Purified Polyclonal Antibody (Catalog # AF3628) (filled histogram) or a Normal Goat IgG isotype control (Catalog # AB-108-C, open histogram) followed by a PE-conjugated Donkey Anti-Goat Secondary Antibody (Catalog # F0107).




**Enrichment of MSCs from C57/BL6 Mouse Compact Bone Using the MagCelect MSC Isolation Kit.** Cells before (A) and after (B) mesenchymal stem cell enrichment were stained with an Alexa 488-conjugated Anti-Mouse PDGF Rα Antibody (clone 189208) and an APC-conjugated Rat Anti-Mouse CD45 Monoclonal Antibody (Catalog # FAB114A).

## MSC Expansion Media

### Features

- **Reliable** – StemXVivo media products are specifically designed and validated for stem cell culture
- **Consistent** – tested for lot-to-lot consistency
- **Flexible** – available as xeno-free, serum-free, or serum-containing media




**Phenotypic Analysis of Human MSCs Expanded in StemXVivo MSC Expansion Media.** Human MSCs were expanded using StemXVivo Mesenchymal Stem Cell Expansion Media (Catalog # CCM004). Filled histograms indicate cells stained with markers of undifferentiated MSCs including PE-conjugated Mouse Anti-Human Endoglin/CD105 Monoclonal Antibody (Catalog # FAB10971P) or PE-conjugated Mouse Anti-Human CD45 Monoclonal Antibody (Catalog # FAB1430). The open histograms show isotype-matched control staining. MSCs appropriately lack expression of CD45.

## CryoDefend-Stem Cells Media

### Features

- **Robust** – greater recovery of viable MSCs compared to conventional media
- **Flexible** – available as serum-free or fully defined, protein free cryopreservation media
- **Validated** – specifically designed and validated for stem cell culture



**Recovery and Marker Expression of Rat MSCs Cryopreserved in CryoDefend-Stem Cells Media.** Rat MSCs (Catalog # PSC003) were cryopreserved in CryoDefend-Stem Cells Media (Catalog # CCM018) or conventional freezing media (95% FBS/5% DMSO) at  $0.7 \times 10^6$  cells/vial. The cryopreserved cells were thawed, cultured for four days, and then assessed for recovery and marker expression. **A)** The number of cryopreserved cells recovered after four days in culture. **B)** CryoDefend preserved MSCs stain positive for APC-conjugated Mouse Anti-Human CD73 Monoclonal Antibody and negative for the PE-conjugated antibodies contained in the Mesenchymal Stem Cell Verification Flow Kit's (Catalog # FMC020) negative marker cocktail.

# Verify


Throughout the expansion and differentiation of MSCs it is important to be confident in the starting populations multipotency. Beginning an experiment with suboptimal, unverified populations will put the investigator at risk for inconsistent results in their downstream experiments, thus wasting time and reagents. R&D Systems offers a series of all-in-one kits that verify MSC multipotency through functional differentiation or cell-specific marker expression.

| Product   | Catalog # |
|---|-----------|
| Human Mesenchymal Stem Cell Functional Identification Kit | SC006     |
| Mouse Mesenchymal Stem Cell Functional Identification Kit | SC010     |
| Rat Mesenchymal Stem Cell Functional Identification Kit   | SC020     |
| Human Mesenchymal Stem Cell 4-Color Flow Kit              | FMC002    |
| Human Mesenchymal Stem Cell Verification Flow Kit         | FMC020    |
| Mouse Mesenchymal Stem Cell 4-Color Flow Kit              | FMC003    |

## Mesenchymal Stem Cell Functional Identification Kits

### Features

- **Reliable** – induces MSC trilineage differentiation with kit-provided supplements
- **Complete** – contains antibodies to confirm successful differentiation
- **Compliant** – defines human MSCs according to International Society for Cellular Therapy (ISCT) recommendations
- **Flexible** – available for verification of human, mouse, and rat MSCs




**Verification of Multipotency using the Human Mesenchymal Stem Cell Functional Identification Kit.** Human MSCs were cultured in StemXVivo Mesenchymal Stem Cell Expansion Media (Catalog # CCM004) and differentiation was induced as indicated using the media supplements included in the Human Mesenchymal Stem Cell Functional Identification Kit (Catalog # SC006). The kit also contains a Goat Anti-Mouse FABP-4 Antigen Affinity-purified Polyclonal Antibody (adipocytes), a Mouse

Anti-Human Osteocalcin Monoclonal Antibody (osteocytes) and a Goat Anti-Human Aggrecan Antigen Affinity-purified Polyclonal Antibody (chondrocytes) for the confirmation of differentiation status. The cells were stained using the NorthernLights™ 557-conjugated Donkey Anti-Goat (Catalog # NL001; red) or Anti-Mouse (Catalog # NL007; red) IgG Secondary Antibodies, and the nuclei were counterstained with DAPI (blue).

## MSC Multi-Color Flow Cytometry Kits

### Features

- **Efficient** – simultaneously detect established positive and negative multipotency markers
- **Flexible** – available for verification of human and mouse MSCs
- **Compliant** – defines human MSCs according to International Society for Cellular Therapy (ISCT) recommendations




**Human Bone-Marrow Derived Cells Fulfill the ISCT's Definition of Human MSCs Based on Marker Expression.** Human bone marrow-derived MSCs were stained using the antibodies and reagents provided in the Human Mesenchymal Stem Cell Marker Verification Flow Kit (Catalog # FMC020). The data shows positive expression of MSC-associated surface antigens CD73, CD90, and CD105. In contrast, minimal expression of antigens recognized by the Negative Marker Cocktail was detected.

## Individual Antibodies for MSC Characterization

Customize your MSC characterization by choosing your preferred R&D Systems Antibodies.

### Our Antibodies:

- **Reliable** – specifically bind to lineage-committed bone marrow-derived cells
- **Flexible** – can be used with magnetic separation systems or with flow cytometry cell sorting for enrichment of uncommitted MSCs



**Detection of Endoglin/CD105 and 5'-Nucleotidase/CD73 on Human MSCs.** Human bone marrow-derived MSCs were stained for positive MSC markers using a APC-conjugated Mouse Anti-Human Endoglin/CD105 Monoclonal Antibody (Catalog # FAB10971A) and a PE-conjugated Mouse Anti-Human 5'-Nucleotidase/CD73 Monoclonal Antibody (Catalog # FAB5795P). Quadrants were set based on isotype controls.

# Differentiate


Efficiently and consistently driving MSCs into adipocytes, chondrocytes, or osteocytes is essential for experimental productivity as well as for reducing costs and labor associated with the lengthy differentiation process. These challenges are remedied by R&D Systems StemXVivo line of base media and supplements, which contain defined, premium quality recombinant proteins to effectively drive MSC differentiation while limiting experimental variation.

## StemXVivo Media and Supplements

### Features

- **Flexible** – available for induction of adipocytes, chondrocytes, or osteocytes
- **User-defined** – combine specific base media and StemXVivo supplements
- **Optimized** – media and supplements developed to support MSC differentiation
- **Diverse** – differentiates MSCs from human, mouse, and rat

| Product                          | Catalog #      |
|----------------------------------|----------------|
| Osteogenic/Adipogenic Base Media | CCM007         |
| Chondrogenic Base Media          | CCM005         |
| Adipogenic Supplement            | CCM011         |
| Chondrogenic Supplements         | CCM006, CCM020 |
| Osteogenic Supplements           | CCM008, CCM009 |



**Rat MSCs were Differentiated into Adipocytes using StemXVivo Adipogenic Media and Supplements.** Rat MSCs were differentiated with StemXVivo Osteogenic/Adipogenic Base Media and StemXVivo Adipogenic Supplement. Adipocyte differentiation was confirmed by staining with a Goat Anti-Mouse FABP4 Antigen Affinity-purified Polyclonal Antibody (Catalog # AF1443) followed by a NorthernLights 557-conjugated Donkey Anti-Goat Secondary Antibody (Catalog # NL001; red). Nuclei were counterstained with DAPI (blue).

## Small Molecules for Differentiation

### Features

- **Precise** – Gain temporal control of differentiation pathways
- **Targeted** – Modulate cell fate by targeting specific signaling pathways
- **Defined** – Minimize the use of animal-derived factors

| Compound                    | Use in Stem Cell Research  | Catalog # |
|-----------------------------|--|-----------|
| Indomethacin                | Induces differentiation of MSCs into adipocytes                      | 1708      |
| Dexamethasone               | Induces differentiation of human MSCs into adipocytes                | 1126      |
| IBMX                        | Induces differentiation of MSCs into adipocytes                      | 2845      |
| Insulin (human) recombinant | Promotes adipogenesis  | 3435      |
| Kartogenin                  | Potently induces chondrogenesis in MSCs                              | 4513      |
| AICAR                       | Induces differentiation of bone marrow-derived MSCs into osteoblasts | 2840      |
| GSA 10                      | Induces differentiation of MSCs into osteoblasts                     | 4918      |
| Purmorphamine               | Induces differentiation of MSCs into osteoblasts                     | 4551      |
| SP 600125                   | Prevents differentiation of MSCs into osteocytes                     | 1496      |
| Zebularine                  | Induces cardiomyogenesis in MSCs                                     | 2293      |
| 5-Azacytidine               | Induces cardiomyogenesis in MSCs                                     | 3842      |
| Pioglitazone hydrochloride  | Improves cardiomyogenesis of MSCs                                    | 4124      |



# Investigate

After carefully validating MSC multipotency, efficiently expanding the cell population, and driving differentiation toward a desired cell lineage it is time to investigate the function of the terminally differentiated cells. Our Proteome Profiler™ Array Kits expedite protein analysis by eliminating the time-consuming gel electrophoresis and protein transfer steps required for a Western blot. Adipogenic, osteogenic, and chondrogenic analytes can be easily quantified using R&D Systems ELISAs, including the Quantikine® and DuoSet Development Kits®.

## Proteome Profiler Antibody Array Kits

### Features



- **Rapid** – analyze the expression level of dozens of cytokines simultaneously
- **Economical** – contains 4 membranes – each cytokine is spotted in duplicate
- **Convenient** – no specialized equipment is required

| Product  | Catalog # |
|--|-----------|
| Proteome Profiler Human Adipokine Array Kit                | ARY024    |
| Proteome Profiler Mouse Adipokine Array Kit                | ARY013    |
| Proteome Profiler Rat Adipokine Array Kit                  | ARY016    |
| Proteome Profiler XL Human Cytokine Array for Osteogenesis | ARY022    |

## Quantikine ELISA and DuoSet Immunoassay Kits

Available for proteins associated with adipogenesis, chondrogenesis, and osteogenesis.

|                        |         | Catalog # |            |
|------------------------|---------|-----------|------------|
|                        | Species | DuoSet    | Quantikine |
| <b>Adipocytes</b>      |         |           |            |
| Adiponectin/Acrp30     | human   | DY1065    | DRP300     |
|                        | mouse   | DY1119    | MRP300     |
|                        | rat     |           | RRP300     |
| Leptin                 | human   | DY398     | DLP00      |
|                        | mouse   | DY498     | MOB00      |
|                        | rat     | DY498     | MOB00      |
| <b>Chondrocytes</b>    |         |           |            |
| Aggrecan               | human   | DY1220    |            |
| SPARC/Osteonectin      | human   |           | DSP00      |
| <b>Osteocytes</b>      |         |           |            |
| Osteocalcin            | human   |           | DSTCNO     |
| Osteopontin            | human   | DY1433    | DOST00     |
|                        | mouse   | DY441     | MOST00     |
|                        | rat     |           | MOST00     |
| Pro-Collagen I alpha 1 | human   | DY622-05  |            |
| SPARC/Osteonectin      | human   |           | DSP00      |



**Expression of Osteoblast-Associated Proteins during Osteogenesis.** Human bone marrow-derived MSCs were differentiated in culture using StemXVivo Osteogenic/Adipogenic Base Media (Catalog # CCM007) and StemXVivo Human Osteogenic Supplement (Catalog # CCM008). Nine days after the induction of osteogenesis, cell culture supernatants were used to assess expression of osteoblast-associated proteins using the Proteome Profiler Human XL Cytokine Array. The data show an expected increase in the expression of KGF/FGF-7, IGFBP-2, GDF-15, CCL5/RANTES, and MIF as cells differentiated into osteoblasts. Additionally, the data demonstrate that Chitinase 3-like 1 and Dkk-1 are strongly expressed in both MSCs (top) and osteoblasts (bottom).

**R&D systems**  
a biotechnique brand

USA & Canada R&D Systems, Inc.  
614 McKinley Place NE  
Minneapolis, MN 55413, USA  
TEL: (800) 343-7475 (612) 379-2956  
FAX: (612) 656-4400  
E-MAIL: [info@RnDSystems.com](mailto:info@RnDSystems.com)  
[RnDSystems.com](http://RnDSystems.com)

R&D Systems Europe Ltd.  
19 Barton Lane, Abingdon Science Park  
Abingdon OX14 3NB, UK  
TEL: +44 (0)1235 529449  
FAX: +44 (0)1235 533420  
E-MAIL: [info@RnDSystems.co.uk](mailto:info@RnDSystems.co.uk)  
[RnDSystems.com](http://RnDSystems.com)

R&D Systems China Co., Ltd.  
24A1 Hua Min Empire Plaza  
726 West Yan An Road, Shanghai, PRC 200050  
TEL: +86 (21) 52380373  
FAX: +86 (21) 52371001  
E-MAIL: [info@RnDSystemsChina.com.cn](mailto:info@RnDSystemsChina.com.cn)  
[RnDSystemsChina.com.cn](http://RnDSystemsChina.com.cn)

International Distributors Please visit [RnDSystems.com/Distributors](http://RnDSystems.com/Distributors) for a full list of international distributors.