




# Biosimilar는 GMP와 함께합니다

## GMP (Good Manufacturing Practice)

- 우수 의약품을 제조하기 위한 공장에 원료의
- 구입~출하까지의 기준 각국마다 자체적인 GMP
- 가이드라인이 있으며 기준을 따르지 않은 경우 불이익 발생

\* 생산단계 뿐만 아니라 연구단계에서도 GMP의 규정은 지켜져야 합니다.



MFDS(식품의약품안전청)에서 제공하는 GMP 월간소식지를 통해 GMP 법령 개정안을 체크해보세요.

[Click here >>](#)

## R&D systems GMP 제품에서 제공하는 제품 정보 항목들



### Recombinant Human TGF-β1 GMP

Catalog Number: 240-GMP

DESCRIPTION	
Source	Chinese Hamster Ovary cell line, CHO-derived human TGF-beta 1 protein Ala279-Ser390 Accession # P01137 Manufactured and tested under cGMP guidelines.
N-terminal Sequence Analysis	Ala-Leu-Asp-Thr-Asn-Tyr-(Cys)-Phe-Ser-Ser
Structure / Form	Disulfide-linked homodimer
Predicted Molecular Mass	12.8 kDa (monomer)
SPECIFICATIONS	
SDS-PAGE	12 kDa, reducing conditions 24 kDa, non-reducing conditions
Activity	Measured by its ability to inhibit the IL-4-dependent proliferation of HT-2 mouse T cells. Tsang, M. <i>et al.</i> (1995) Cytokine 7:389. The ED <sub>50</sub> for this effect is 0.04-0.2 ng/mL. The specific activity of recombinant human TGF-β1 GMP is approximately 2.5 x 10 <sup>4</sup> U/μg, which is calibrated against human TGF-β1 Standard (NIBSC code: 89/514).
Endotoxin Level	<0.10 EU per 1 μg of the protein by the LAL method.
Purity	>97%, by SDS-PAGE with silver staining.
Host Cell Protein	< 0.5 ng per μg of protein when tested by ELISA.
Mycoplasma	Negative when tested in a ribosomal RNA hybridization assay.
Adventitious Virus	Master Cell Bank tested for adventitious viruses
Formulation	Lyophilized from a 0.2 μm filtered solution in Acetonitrile and TFA. See Certificate of Analysis for details.
PREPARATION AND STORAGE	
Reconstitution	Reconstitute at 100 μg/mL in 4 mM HCl.
Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> <li>• A minimum of 12 months when stored at ≤ -20 °C as supplied. Refer to lot specific COA for the Use by Date.</li> <li>• 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>• 3 months, ≤ -20 °C under sterile conditions after reconstitution.</li> </ul>

- Source information
- N-terminal sequencing of the first 10 amino acids
- Purity by SDS-PAGE and Mass Spectrometry
- Bioactivity testing calibrated to WHO International Standards
- Sterility testing to USP
- Endotoxin levels
- Certification and Regulatory Guidelines followed, including:
  - ▶ ISO 9001:2015, ISO 13485:2016-certified facility
  - ▶ USP Chapter <1043>, Ancillary Materials for Cell, Gene, and Tissue-Engineered Products
  - ▶ USP Chapter <92>, Growth Factors and Cytokines Used in Cell Therapy Manufacturing
  - ▶ Ph. Eur. General Chapter 5.2.12, Raw Materials of Biological Origin for the Production of Cell-based and Gene Therapy Medicinal Products