
Toxicology - CometAssay

R&D SYSTEMS
a **biotechne** brand

Comet Assay Procedure

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In vitro assessment of the photo(geno)toxicity associated with Lapatinib, a Tyrosine Kinase inhibitor

[Click here >>](#)


The epidermal growth factor receptors EGFR and HER2 are the main targets for tyrosine kinase inhibitors (TKIs). The quinazoline derivative lapatinib (LAP) is used since 2007 as dual TKI in the treatment of metastatic breast cancer and currently, it is used as an oral anticancer drug for the treatment of solid tumors such as breast and lung cancer. Although hepatotoxicity is its main side effect, it makes sense to investigate the ability of LAP to induce photosensitivity reactions bearing in mind that BRAF (serine/threonine-protein kinase B-Raf) inhibitors display a considerable phototoxic potential and that afloqualone, a quinazoline-marketed drug, causes photodermatitis. Metabolic bioactivation of LAP by CYP3A4 and CYP3A5 leads to chemically reactive N-dealkylated (N-LAP) and O-dealkylated (O-LAP) derivatives. In this context, the aim of the present work is to explore whether LAP and its N- and O-dealkylated metabolites can induce photosensitivity disorders by evaluating their photo(geno)toxicity through in vitro studies, including cell viability as well as photosensitized protein and DNA damage. As a matter of fact, our work has demonstrated that not only LAP, but also its metabolite N-LAP have a clear photosensitizing potential. They are both phototoxic and photogenotoxic to cells, as revealed by the 3T3 NRU assay and the comet assay, respectively. By contrast, the O-LAP does not display relevant photobiological properties. Remarkably, the parent drug LAP shows the highest activity in membrane phototoxicity and protein oxidation, whereas N-LAP is associated with the highest photogenotoxicity, through oxidation of purine bases, as revealed by detection of 8-Oxo-dG.

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STEP1

Electrophoresis system

- Electrophoresis tank.
- power supply.

+


STEP2

CometSlide, Cometchip

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- 2well
- 20 well
- 96 well

↓

STEP3

Data analysis

전용 Automated software

- Cometanalysis program