

Upgrade Your **Western Blot**

Gel-free, Blot-free, Hands-free, **Auto-Western blot**



High Sensitivity: 0.38pg까지 Detection 가능

Less Sample & less Ab: 0.6ug Sample (최소농도 0.25mg/ml)

Multiplex & Ab Striping: 최대 8개 targets 동시 분석 / Sample

3hrs Run time: from loading to result

Total Protein Normalization: 어떤 실험 조건에서도 흔들리지 않는 Normalization

Full Automation: 높은 재현성

1,000개 이상 Certified BioTechne (R&D + Novus) Abs

웅비 Assay Service 의뢰 가능


장비 데모 가능 : 웅비실험실(1일)




Wes

Features & Advantages


01 High-Sensitivity



02 Total Protein Normalization



03 Ab Striping & Multiplex



04 데모 or 분석서비스 가능



05 필요한 스펙에 맞춰 선택하세요!



RUN TIME	3 HOURS	
SAMPLES	12-24	
SEPARATION	SIZE	
DETECTION	CHEMILUMINESCENCE FLUORESCENCE	CHEMILUMINESCENCE
RePlex™ (Strip & ReProbe Replacement)	YES	NO
OTHER CAPABILITIES	BLOT IMAGING	

응비에서 분석서비스 받고
편안하게 연구하세요!

