


+ The Best

Cytotoxicity LDH Assay Kit




- ✓ **Damaged Cell**에서 분비되는 LDH 측정을 통한 **Cytotoxicity** 분석
(Cell toxicity관련, MTT Data 보강 /보완에 최적)
- ✓ **Homogeneous & Non-Homogeneous** 분석
(attached cells) (suspended cells)
- ✓ **Longer Working Solution Stability**
(Storage Condition : 0-5 °C for 6 months)
- ✓ **Ready to use**

왜 LDH Kit를 사용 할까요?

Measurement of Live Cells


Live cells: metabolic activity(CCK-8, MTT, MTS, XTT etc)



Measurement of Live and Dead Cells


Live cells: metabolic activity(CCK-8, MTT, MTS, XTT etc)

Dead cells: destruction of plasma membrane(LDH Assay Kit)



누가 LDH Kit를 사용 할까요?

Dojindo LDH를 구입한 새 고객 중 CCK-8고객의 퍼센티지(2016/6-2017/5)



세포 독성 검사를 위한 LDH Assay Kit

