

## Quickkit™ ELISAs

Get Quantikine Quality in



\*출시기념 이벤트: QuickKit 전품목 **499,000원** /1 plate (2020.5/1 ~ 7/31)

### ✓ FAST

> 90분 안에 결과 확인 할 수 있는 가장 빠른 ELISA





### ✓ Only 1 Wash Step

> 간단하고, 편리하게



### ✓ 정확 & 재현성 있는 결과

> R&D Quantikine만의 diluents 조합






	QuicKit ELISA 종류	Range (pg/mL)	행사가격	Cat. No.
Human	Total IL-18	15.6 - 1,000	499,000	QK318
	CRP	800 - 50,000	499,000	QK1707
	IL-2	31.3 - 2,000	499,000	QK202
	IL-6	12.5 - 800	499,000	QK206
	TNF-alpha	31.3 - 2,000	499,000	QK210
	G-CSF	39.1 - 2,500	499,000	QK214
	TNF RI/TNFRSF1A	15.6 - 1,000	499,000	QK225
	CXCL10/IP-10	15.6 - 1,000	499,000	QK266
	CD14	188 - 12,000	499,000	QK383
	Leptin	15.6 - 1,000	499,000	QK398

## Assay Service

Luminx / ELISA / Western  
Ab Array / Ella



High Quality Time Saving!  
웅비메디텍 Assay Service를 통해  
합리적인 가격으로 신속 정확한 서비스를 만나세요

## Table Seminar

지금 웅비메디텍에 신청하세요



- Tool for Cytotoxicity
- Immunoassay Selection Guide
- 3D Cell culture & Analysis Tools
- New Technologies for Western blot
- ELISA Technical Tip & QnA
- in vitro Drug Screening Tools
- DNA Damage & Apoptosis