

IHC 연구자들을 위한 빠르고 선명한

VisUCyte HRP Polymer Antibody

HRP-polymer based signal amplification

VisUCyte만의 장점

- ✔ Sensitivity 30-50배 증가
- ✔ Only 2 step
- ✔ Primary Ab 사용량 최대 1/3 감소
- ✔ Background 최소화
- ✔ 비용 절감 & 시간 단축



Order Information

	#Cat
Goat anti-Mouse IgG	VC001
Goat anti-Mouse/Rabbit IgG	VC002
Goat anti-Rabbit IgG	VC003
Donkey anti-Goat IgG	VC004
Goat anti-Rat IgG NEW	VC005
Donkey anti-Sheep IgG NEW	VC006




Related Product

	#Cat
Antigen Retrieval Reagent-Basic	VCTS021
Antigen Retrieval Reagent-Acidic	VCTS022
Antigen Retrieval Reagent-Universal	VCTS023

Polymer-Based VS ABC 비교 Data


Polymer backbone
 Avidin/Streptavidin
 Biotin
 HRP
 Secondary antibody
 Primary antibody
 Antigen

Amplification





Polymer-Based

Amplification



ABC

+ primary antibody



- primary antibody

