



항체신약 & Biosimilar 개발의 필수 Tools Fc gamma Receptor Proteins



IgG Binding ability check 완료!




Coating: FCRN (8639-FC)
Detection: Biotinylated human IgG
ED₅₀: 2-10 µg/mL.



Immobilized: Fc gamma RIIIA/CD16a (4325-FC)
Sample: HumanIgG
K_d < 50nM

SDS-PAGE를 이용한 purity check 완료!











Fc gamma RIIA/CD32a Avi-tag His-tag Protein (AVI1330) Purity >95%

Cat#	Name	Isoform	Tag
1257-FC	Fc gamma RI/CD64		His
4325-FC	Fc gamma RIIA/CD16a	V176	His
9595-CD	Fc gamma RIIA/CD32a	H167	His
1330-CD	FC gamma RIIA/CD32a	R167	His
AVI1330	Fc gamma RIIA/CD32a		Avi, His
1875-CD	Fc gamma RIIB/C (CD32b/c)		His
AVI1875	Fc gamma RIIB/CD32b		Avi
AVI8894	Fc gamma RIIIA	F176	His, Avi
8894-FC	Fc gamma RIIIA/CD16a	V176F	His
AVI10468	Fc gamma RIIIA/CD16a		His, Avi
1597-FC	Fc gamma RIIIB/CD16b		His
9948-FC	Fc gamma RIIIB/CD16b (aa18-200)		His
8639-FC	FCRN		His

Anti-Idiotypic Antibody

- Binding assay
- Flow Cytometry
- Binding Inhibition
- Direct ELISA



Anti-Idiotypic Antibody	Uncojugated Abs	Available Flow Cytometry Fluorochrome
Anti- Cetuximab	MAB9626 (ELISA, Flow Cytometry, Simple Western, Western Blot)	Alexa Fluor 350, 405, 488, 594, 647, 700, 750
Anti- Rituximab	MAB9630 (Block/Neutralize, ELISA, Flow Cytometry, Western Blot)	Alexa Fluor 350, 405, 488, 594, 647, 700, 750
Anti- Trastuzumab	MAB95471 (ELISA, Flow Cytometry, Western Blot) MAB9547 (ELISA, Western Blot)	Alexa Fluor 350, 405, 488, 594, 647, 700, 750
Anti- Adalimumab	MAB9546 (ELISA, Western Blot) MAB9616 (ELISA, Western Blot)	
Anti- Alemtuzumab	MAB10246 (Flow Cytometry)	Alexa Fluor 350, 405, 488, 594, 647, 700, 750
Anti- Gemtuzumab	MAB96171 (ELISA, Flow Cytometry, Western Blot) MAB9617 (ELISA, Western Blot)	Alexa Fluor 350, 405, 488, 594, 647, 700, 750
Anti- Aducanumab	MAB9545 (ELISA, Western Blot) MAB10180 (Western Blot)	
Anti- Teprotumumab	MAB96341 (Block/Neutralize, Flow Cytometry) MAB9634 (ELISA, Western Blot) MAB96342 (Western Blot)	Alexa Fluor 350, 405, 488, 594, 647, 700, 750
Anti- Teplizumab	MAB10245 (Flow Cytometry, Western Blot)	Alexa Fluor 350, 405, 488, 594, 647, 700, 750
Anti- Brentuximab	MAB9584 (ELISA, Flow Cytometry, Western Blot)	Alexa Fluor 350, 405, 488, 594, 647, 700, 750
Anti- Solanezumab	MAB10124 (CyTOF-ready, Flow Cytometry, Western Blot)	
Anti- Daclizumab	MAB10218 (Western Blot)	
Anti- Ibritumomab	MAB9544 (ELISA, Western Blot)	
Anti- Ipilimumab	MAB10064 (Western Blot)	
Anti- Satumomab	MAB10219 (Western Blot)	