

NSC 연구자들을 위한 Culture Total Solution



Culture Supplement

B27 대비 더 효율적인 cell 분화 검증!

효율적인 cell 분화를 원하신다면?

R&D의 Supplement

Good Quality | Low Toxicity | Efficiency

Product	Catalog #
N-2 Plus Media Supplement	AR003
N-2 MAX Media Supplement	AR009
N21-MAX Media Supplement	AR008
N21-MAX Insulin Free Media Supplement	AR010
N21-MAX Vitamin A Free Media Supplement	AR012



Small Molecules

믿을 수 있는 Quality! Tocris Small Molecules!

High Purity

다양한 제품군

믿을 수 있는 QC 분석

저렴한 가격

Product	Catalog #	Product	Catalog #
17-AAG	1515	INDY	4997
Amiodarone hydrochloride	4095	P7C3	4076
CHIR 99021	4423	SB 431542	1614
Cyclopamine	1623	SB 203580	1202
DMH-1	4126	Y-27632 dihydrochloride	1254



Growth Factor

Quality 좋기로 소문난 R&D Growth Factor로 더 만족스러운 결과를 얻으세요!

Recombinant Protein	Catalog #			ProDots Proteins, Catalog #	GMP Grade	Recombinant Protein	Catalog #			ProDots Proteins, Catalog #	GMP Grade
	Human	Mouse	Rat				Human	Mouse	Rat		
BDNF	248-BD	248-BD	248-BD		248-GMP	IGF-I	291-G1	791-MG	4326-RG	PRD291	291-GMP
β -NGF	256-GF	256-GF	556-NG			Noggin	6057-NG	1967-NG			6057-GMP
BMP-2	355-BM	355-BM	355-BM		355-GMP	NT-4	268-N4	3236-N4			268-GMP
BMP-4	314-BP	5020-BP		PRD314	314-GMP	PDGF	120-HD				
CNTF	257-NT		557-NT			PDGF-AA	221-AA		1055-AA		
EGF	236-EG	2028-EG	3214-EG	PRD236	236-GMP	PDGF-AB	222-AB		1115-AB		
FGF basic	233-FB	3139-FB	3339-FB	PRD233	233-GMP	PDGF-BB	220-BB		520-BB	PRD220	220-GMP
FGF-4	235-F4	5846-F4				Sonic Hedgehog/Shh		464-SH			1314-GMP
FGF-8b	423-F8	423-F8				VEGF 165	293-VE			PRD293	293-GMP
GDNF	212-GD		512-GF								


Rat Cortical Neuron Culture Kit

- High Quality R&D BDNF, IGF-1 제공
- Rat cortical neuron culture에 최적화 된 Kit
- Embryonic / Postnatal 모두 가능

COMPONENTS	PART #	# VIALS
Neuronal Base Media	967977	1 vial (500 mL)
N21-MAX Media Supplement (50X)	893414	1 vial (10 mL)
Recombinant human BDNF	967978	1 vial
Recombinant human IGF-I	967979	1 vial
Reconstitution Buffer 1	967552	2 vials (1 mL/vial)



R&D competitor



Good Quality & Efficiency

똑같은 culture time, 더 많은 neurite outgrowth와 synaptic puncta! R&D Kit의 Quality와 Efficiency를 보여주는 결과입니다!

* Provided this within the expiration date of the kit.