

Cyto3D™

Live-Dead Assay Kit




3D Live-Dead Assay에 최적화된 KIT!

- Ready-to-use
- One-Pipetting으로 **Two-Targeting**
- 3D 뿐만 아니라 2D도 사용 가능!
- 5- 10분의 소요시간
- 냉장에서 15개월까지 안정적!



Easy Setup and Use




1 Add Cyto3D

2 Ready for detection

Specifications

Formulation	Mixed Acridine Orange(AO) and PI
Detection Method	Fluorescent
Excitation/Emission	AO(494/517nm), PI(535/617nm)
Storage	2 ~ 8 °C



Matrix내의 Spheroids를 Live Image로 확인 가능!

Order Information

Product	Cat#	Size
Cyto3D™ Live-Dead Assay Kit	BM01	1 mL (500 test (at 2 μL per 100 μL)

