


VitroGel HEK293






HEK293 Cell Scale Up Culture !

- CHO Cell 말고 human cell HEK293을 활용한 단백질 생산
- Stirring culture가 아닌 3D Static culture 로 생산 UP!

- Xeno-free
- Ready to Use
- Easy Cell Harvest
- 자동화 장비 이용 가능
- 실온에 안정
- Injectable
- Biosafety
- Mix & Match



▶ 3D static suspension culture



A. 3D static suspension culture of HEK293 cells in VitroGel HEK293 The cell proliferation was tested by Cell Counting Kit-8 (CCK-8) assay.

B. Cell suspension can be prepared at high cell density (10^7 cells/mL) to reach the peak of protein expression within 48 to 72 hours.

