

# 15주년 기념 핑큐 세일

# SALE

15.08.24(월)~15.09.30(수)

Dojindo 베스트 아이템  
Cell Counting Kit-8

## 단독 30% 할인

- ◆ 1,000 T (CK04-11 : 5 ml x 2 bottles)
- ◆ 3,000 T (CK04-13 : 5 ml x 6 bottles)
- ◆ 10,000 T (CK04-20 : 50 ml x 2 bottles)

놓치지 않을 거예요~



제품 주문 및 상담: [techserv@woongbee.com](mailto:techserv@woongbee.com)/031-776-3300



# The Best for Cell Proliferation & Cytotoxicity

# CCK-8



- ◆ More sensitive than MTT, XTT, MTS, WST-1
- ◆ No Toxicity to Cells
- ◆ Ready-to-use/ one-bottle
- ◆ Long Stability (냉장 1년)
- ◆ Less sensitive to light
- ◆ No radioisotope or DMSO
- ◆ No harvest, washing or solubilization

MTT...  
최선일까요?




CCK-8 vs. MTT 비교동영상




## Most simple procedure

Procedure	Preparation	Add	Dissolve	Measure
CCK-8		✓		✓
MTS WST-1	✓	✓		✓
MTT	✓	✓	✓	✓

## No Toxicity to Cells



CCK-8




MTS




MTT


## Highest Sensitivity



## Longer stability



## Highest Cited Kit



신제품

# LDH Assay for Drug toxicity

Drug 처리 후에 손상되거나 죽은 세포로부터 Drug Toxicity 측정

- ✓ Homogeneous, non-homogeneous assays 모두 가능
- ✓ Stable (0~5 degree) - not need to freeze
- ✓ 사용하기 편리한 colorimetric assay




Fig.1 Principle of cytotoxicity measurement