

Notch Signaling 관련 TOCRIS Products



## Products for Notch Signaling

Cat.#	Name	Description
<b>Modulators</b>		
4898	TC-E 5006	y-secretase modulator; reduces A $\beta$ 42 levels
<b>Inhibitors</b>		
5172	FLI 06	Inhibitor of Notch Signaling
4283	Begacestat	y-secretase inhibitor; lowers A $\beta$ 40 and A $\beta$ 42 levels
2870	BMS 299897	Potent y-secretase inhibitor
2654	Compound W	y-secretase inhibitor
2634	DAPT	y-secretase inhibitor
4489	DBZ	y-secretase inhibitor; inhibits Notch pathway
4495	Flurizan	y-secretase inhibitor; lowers A $\beta$ 42 levels in vitro
2677	JLK 6	Inhibitor of y-secretase-mediated $\beta$ APP processing
2627	L-685,458	Potent and selective y-secretase inhibitor
4000	MRK 560	y-secretase inhibitor; attenuates amyloid plaque deposition

