



Products for AKT Signaling



PI3-Kinases

Cat.#	Name	Description	Cat.#	Name	Description
Inhibitors					
5595	A66	Potent and selective PI3-kinase p 110α inhibitor	2418	LY 303511	Negative control of LY 294002(Cat. No.1130)
3671	AS 252424	Selective inhibitor of PI3-kinase γ	3977	3-Methyladenine	Class III PI 3-kinase inhibitor, also inhibits autophagy
3578	AS 605240	Potent and selective PI3-kinase γ(PI3-Kγ)inhibitor	4820	PF 04691502	Potent and selective dual PI 3-K/mTOR inhibitor
4839	AZD 6482	Potent and selective PI3-Kβ inhibitor	4823	PF 05212384	Potent and selective dual PI3-K/mTOR inhibitor
3606	BAG 956	Dual PI 3-kinase and PDK1 inhibitor	2930	PI 103 hydrochloride	Inhibitor of PI 3-kinase, mTOR and DNA-PK
4674	CZC 24832	Selective inhibitor of PI 3-kinase γ	2814	PI 828	PI 3-kinase inhibitor, more potent than LY 294002 (Cat. No. 1130)
4702	ETP 45658	Potent PI3-kinase inhibitor; also inhibits DNA and mTOR	3894	PP 121	Dual kinase inhibitor; inhibits PI 3-K family kinases
4026	GSK 1059615	Potent PI3-kinase inhibitor	1125	Quercetin	Non-selective PI 3-kinase inhibitor
4840	KU 0060648	Dual DNA-PK and PI3-K inhibitor	4264	TG 100713	PI 3-kinase inhibitor
1130	LY 294002 hydrochloride	Prototypical PI 3-kinase inhibitor, also inhibits other kinases	1232	Wortmannin	Potent, irreversible inhibitor of PI 3-kinase. Also inhibitor of PLK1
Activators			Other		
1983	740 Y-P	Cell-permeable PI3-kinase activator	5126	iMDK	Suppressor of MDK expression; inhibits PI 3-K signaling