

# illumigene™ *C. difficile*

Molecular Simplified

**A New Standard of patient care**  
**Definitive Answers, Confident Results**

non-PCR (based on **LAMP**)

The **Most Accurate**

Toxigenic Culture	<b>Meridian illumigene</b>	Cepheid Xpert	VIDAS Toxin A & B	Toxin A & B EIA
Sensitivity	<b>95.2%</b>	93.5%	80.0%	80.8%
Specificity	<b>95.3%</b>	94.0%	97.5%	97.3%
PPV	<b>92.5%</b>	74.0%	76.4%	78.0%
NPV	<b>99.5%</b>	98.8%	97.8%	97.9%

FDA/ NHS 자료 참조

Save **Healthcare costs**

No DNA Extraction

TAT < **40 min**



FDA Cleared

**Group A Strep., B Strep.**  
***Mycoplasma pneumoniae*\***

