

# Introducing Tocriscreen Plus™

A unique new library of 1,280 biologically active compounds


Tocris® has over 10 years' experience of producing compound libraries and this, combined with feedback from our customers and the latest advances in research, adds up to a unique new library, the Tocriscreen Plus™ / Tocriscreen Plus Mini™

## Tocriscreen Key Highlights

- ▶ HTS와 HCS에 적합한 biologically active compound collections
- ▶ 1280개의 compounds 구성
- ▶ 98% 이상의 고순도 compounds
- ▶ Ready to use 10mM DMSO
- ▶ Matrix tube와 SeptraSeal caps를 포함한 96well racks 형태
- ▶ Custom Compound Library Service 가능
- ▶ **NEW!** FDA-Approved Compounds 이용 가능 (160개)



**Tocriscreen Plus Targets**



**TOCRISCREEN™**  
PLUS

**TOCRISCREEN™**  
PLUS Mini

**TOCRISCREEN™**  
FDA-Approved Compounds

Category	Cat. No.	Product Name	Description	Unit Size
<b>Tocriscreen Compound Libraries</b>				
	5840	Tocriscreen Plus	A library of 1280 biologically active compounds (250µL 10mM DMSO solutions)	1 Set
	5841	Tocriscreen Plus Mini	A library of 1280 biologically active compounds (50µL 10mM DMSO solutions)	1 Set
<b>NEW!</b>	5932	Tocriscreen FDA-Approved	A library of 160 biologically active compounds (250µL 10mM DMSO solutions)	1 Set