

60분만에 COVID-19-CRS 확인

Automated high Sensitive Multiplex ELISA

Ella

Extremely Simple: 3 step & 1시간 내 분석완료

High Sensitivity & Wide Dynamic Range: 1 ~ 10,000 pg/mL

비교불가 재현성: No Standard curve (Factory Calibrated)


High Correlation w/Quantikine ELISA kit

25 uL Sample vol. 으로 Triplicate 분석

다양한 cartridge:

- Single Plex
- Multiplex
- Customizable

웅비 Assay Service 의뢰가능 ☆



COVID-19 Cytokine Analysis Panel

Plex#	Analytes
4-plex	IL-1 β , IL-6, IL-8/CXCL8, TNF- α
4-plex	IL-7, IL-10, IL-12p70, TNF- α
8-plex	IFN- γ , IL-1 β , IL-2, IL-4, IL-6, IL-10, IL-12p70, TNF- α
8-plex	CCL2/MCP-1, CCL3/MIP-1 α , CXCL10/IP-10, IL-1 β , IL-2, IL-7, IL-10, TNF- α

for rapid multicytokine monitoring of COVID-19 patients
in Mt Sinai Hospital (New York)

연구 자세히 보러가기 ▶▶

