


# biotechne®

## ACE-2 ELISA Kit

**R&D SYSTEMS**  
a biotechne® brand

### Human ACE-2 DuoSet ELISA




- Catalog number : DY933-05
- Standard range : 0.313 – 20 ng/mL
- Sample type : Cell culture supernates
- Ancillary Kit : DY008

[▶ R&D Systems Human ACE-2 DuoSet ELISA 정보 확인하기](#)

**R&D SYSTEMS**  
a biotechne® brand

### Human ACE-2 ELISA



- Catalog number : NBP2-78734
- Standard range : 0.39 – 25 ng/mL
- Sample type : Serum, Plasma,  
and other biological fluids

[▶ Novus Biological Human ACE-2 ELISA 정보 확인하기](#)