




MitoBright (Mitochondrial Staining)

- Mitochondrial Staining
- High Intracellular Retention
- Long-Term Visualization of Mitochondria
- Provided as Ready-to-Use DMSO Solution
- Allows for Collagen-coated Plate


Technical info

	MitoBright LT			Company (T)		
	Green	Red	Deep R	Green	Red	Deep R
Product condition	DMSO solution			Solid		
Intracellular Retention	○	○	○	×	×	×
Stain with serum-containg medium	○	○	○	△	△	△
Fixation after staining(PFA)	○	○	○	○	○	○
Fixation after staining(MeOH)	×	×	×	×	×	×


Fluorescence Properties



High Intracellular Retentivity



Stained in serum-contained media



MitoBright 제품 더보기