

QantiFERON®-TB Gold In-Tube 검사 전 준비사항

장비 및 비품	필수 사항	선택 사항
<u>Blood collection</u>		
<u>Whole blood simulation</u>		
1. 37°C incubator (humidified and CO ₂ not required)	O	
2. 1~1.5ml microtubes 또는 uncoated 96-well plates	O	
3. Biohazard Cabinet Class II		O
<u>ELISA</u>		
1. Distilled water (ELISA grade) (1L)	O	
2. Calibrated pipettes (1ml)	O	
3. 1ml vol. blue tips	O	
4. Calibrated pipettes (~200ul)	O	
5. ~200ul vol. yellow tips	O	
6. Multichannel pipet (50ul,100ul)		O
7. 1~1.5ml microtube 또는 Test tube	O	
8. Reagent reservoirs (polypropylene only)		O
9. Measuring cylinder (1~2L 용량)	O	
10. Microplate shaker	O	
11. Microplate Reader (450nm filter, 620-650nm reference filter)	O	
12. Automated Microplate washer		O
13. Vortex		O
14. Timer		O